



# TOP LTE Quick Release™


## \*\* Therminator™ Chiller Mounting Bracket \*\*

\*\*\*\* PLEASE READ THOROUGHLY PRIOR TO USE FOR IMPORTANT SAFETY INFORMATION \*\*\*\*

### Unpacking:




Please ensure that you have the following items packaged with your mounting bracket.

- 1 Mounting plate
- 1 Mounting bracket
- 2 1/4" bolt
- 2 1/4" nut



### Assembly:

Assembly is quite simple. Place the 1/4" bolts through the mounting plate and LTE stand. Secure with the 1/4" nut (Fig. 1). Mount the chiller to the mounting bracket using the mounting studs on the chiller and the nuts included with the chiller as shown in Fig. 2. Set the mounting bracket into the mounting plate slots as shown in Fig. 3.



### Warning:

Do NOT soak the mounting bracket or mounting plate in sanitizer or cleaning chemicals as this will potentially damage the zinc plating on the bracket. Mild detergent is sufficient to remove soils after your brew day.

©Blichmann Engineering, LLC 2015  
For replacement parts visit [blichmannengineering.com/genuine-replacement-parts](http://blichmannengineering.com/genuine-replacement-parts)

This product uses food grade materials anywhere the product touches the beverage.  
Warning: This product contains or may contain chemical(s) known to the State of California to cause cancer, birth defects, or other reproductive harm.