



# TopTier™


## \*\* Quick Release Therminator™ Bracket \*\*

\*\*\*\* PLEASE READ THOROUGHLY PRIOR TO USE FOR IMPORTANT SAFETY INFORMATION \*\*\*\*

### Unpacking:



Please ensure that you have the following items packaged with your mounting bracket.

- 1 Mounting plate
- 1 Mounting bracket
- 2 5/16" bolt
- 2 5/16" square nut



### Assembly:

Assembly is quite simple. Place the 5/16" bolts through the mounting plate with the 5/16" square nuts. It is helpful to thread the square nuts only a turn or two onto the bolts to make insertion into the t-slot on the stand easier. Slide the bracket assembly into the t-slot and tighten the nut snugly (20 ft-lb or 27 N-m) using a torque wrench (Fig. 1). Caution: excessive tightening may permanently deform the mast. Then mount the chiller to the mounting bracket using the mounting studs on the chiller and the nuts included with the chiller as shown in Fig. 2. Set the mounting bracket into the mounting plate slots as shown in Fig. 3.



### Warning:

Do NOT soak the mounting bracket or mounting plate in sanitizer or cleaning chemicals as this will potentially damage the zinc plating on the bracket. Mild detergent is sufficient to remove soils after your brew day.